

# LC Logger

Low Cost single use  
temperature data logger



## PHARMACEUTICAL

## FRESH

## FROZEN

## TECHNICAL



18 Months  
Shelf Life

Automatic  
PDF & CSV  
Report on  
board

Up to 100  
days recording

Comes with a  
Certificate  
of Validation

Secure and  
easy download  
of encrypted  
data files



### Built-in USB Connector

Directly connect to  
computer,  
no cable required



### IP67 Protection class

Food-grade  
and water-proof  
packaging



### Automatic PDF generation

No software  
required




### Qualified for air transport use

Complies to RTCA/  
DO-160G Environmental  
Conditions, qualified  
for air transport under  
active status



## TEMPERATURE RECORDER



## TECHNICAL SPECIFICATIONS

Device type:	single use internal sensors
Temperature unit:	C or °F
Operating range :	-30C to +70°C
Temperature accuracy	0.5°C
Resolution :	0.1°C/F
Time reference :	UTC+00
Storage capacity :	14, 400 memory points
Recording interval :	Default 10 mn- Factory adjustable
Standard start delay :	Default 30 mn- factory adjustable
Battery:	CR2032
Start/Stop mode :	Push Start/Stop button - Automatic stop when connected to computer
Alarms :	2 alarms in option - Factory adjustable
LED:	RED, GREEN/BLUE-START/STOP-indication
Reporting :	Automatic generation of PDF / Dual scale C F
Computer interface :	USB 2.0
Size :	60 x 20 x 5,7 mm.
Weight :	8g
Alarms :	single event, up to 2 alarms
Certificatlons :	CE, ROHS, UN 38.3
IP rating :	IP67
Validation	with NIST traceable certificate for each unit

## EASY TO MANAGE

**1**

Click the START button to turn on and enable BioXTEND Pro. The data logger is preconfigured and ready to use.

**2**

Can you adapt with your device ? If not, just remove this part

**3**

Remove the device from the delivery and click STOP to disable it. Connect BioXTEND directly to the computer USB port.

**4**

The report (in PDF & CSV format) will download automatically without any software installation.

## BioXTEND Inc.